

ONYX ONE CLINICAL SUMMARY



BACKGROUND

- > ONYX ONE OVERVIEW
- > SHORT-DAPT TRIALS COMPARISON
- > WHY RESOLUTE ONYX™ DES FOR HBR PATIENTS ON 1-MONTH DAPT?



CLINICAL DATA

- > ONYX ONE GLOBAL TRIAL 1-YEAR RESULTS
- > ONYX ONE GLOBAL TRIAL 2-YEAR RESULTS
- > ONYX ONE CLEAR ANALYSIS
- > ACS SUBANALYSIS
- > AF SUBANALYSIS
- > COMPLEX PCI SUBANALYSIS
- > SEX SUBANALYSIS
- > DIABETIC SUBANALYSIS

ONYX ONE OVERVIEW

ONYX ONE GLOBAL TRIAL

First, randomised, 1-month DAPT trial comparing a DES to a DES in high bleeding risk (HBR) patients and proving Resolute Onyx DES safety and efficacy.¹

ONYX ONE CLEAR ANALYSIS

Evaluated Resolute Onyx DES in ~1500 complex HBR patients on 1-month DAPT and reinforced safety and efficacy results from the Onyx ONE Global trial.²

FIRST DES CE INDICATED FOR 1-MO DAPT IN HBR PATIENTS



FIRST DES FDA APPROVED FOR HBR PATIENTS WITH 1-MO DAPT



Real-world, HBR patient population led to **meaningful subanalyses** for greater confidence making short-DAPT decisions

ACS SUBGROUP³
ESC 2020

AF SUBGROUP⁴
TCT 2020

COMPLEX PCI SUBGROUP⁵
TCT 2020

SEX SUBGROUP⁶
SCAI 2021

DIABETIC SUBGROUP⁷
PCR 2021

¹Windecker S, et al. *N Engl J Med*. 2020;382:1208-1218.

²Kandzari D, et al. *Circ Cardiovasc Interv*. 2020;13: e009565.

³Kedhi A, et al. Outcomes in HBR with ACS with 1-Month DAPT. Presented at ESC2020.

⁴Pasupati S, et al. Ischemic and Bleeding Outcomes in Patients With vs. Without AF. Presented at TCT2020.

⁵Kandzari D, et al. Complex PCI with 1-month DAPT in HBR Patients: Presented at TCT2020.

⁶Mehran, R et al. Sex-based Outcomes after PCI in Complex High Bleeding Risk Patients: Results from the Onyx ONE Clear Trial. Presented at SCAI 2021.


⁷Kedhi E, et al. Diabetic High Bleeding Risk Patients with One-Month DAPT: Onyx ONE Clear Results. Presented at Euro PCR 2021.



ONYX ONE GLOBAL TRIAL 1-YEAR RESULTS

Showed Resolute Onyx™ DES was safe and effective in HBR patients on 1-month DAPT

LANDMARK ANALYSIS AFTER DAPT DISCONTINUATION† SHOWED LOW EVENT RATES



- Randomised controlled trial
- Compared Resolute Onyx DES to BioFreedom DCS in ~2000 HBR patients on 1-month DAPT
- 1.6 high bleeding risk criteria per patient
- BARC 3–5 bleeding rate:
 - Resolute Onyx 4.9%
 - BioFreedom DCS 4.4% $p = 0.67$
- Primary endpoint met with Resolute Onyx DES (17.1%) noninferior to BioFreedom DCS (16.9%) for cardiac death (CD), myocardial infarction (MI), and stent thrombosis (ST)
- Resolute Onyx is the first DES CE indicated for 1-month DAPT in HBR patients based on the results from this study

*Third-party brands are trademarks of their respective owners.
 †From 1 month to 1 year. Rates are taken from Kaplan-Meier estimates.
 **Post-hoc analyses were not powered.
 Source: Windecker S, et al. *NEngl J Med*. 2020;382:1208-1218.




Read the Onyx ONE global trial article printed in *The New England Journal of Medicine*




ONYX ONE GLOBAL TRIAL 2-YEAR RESULTS

Reinforced Resolute Onyx™ DES safety and efficacy in HBR patients on 1-month DAPT at two years

SIGNIFICANTLY LOWER SPONTANEOUS MI WITH RESOLUTE ONYX DES†



LOW STENT THROMBOSIS CONFIRMED AT 2-YEAR FOLLOW-UP†



- Final follow-up of the Onyx ONE Global trial at two years
- 1.6 high bleeding risk criteria per patient
- No difference between Resolute Onyx DES (21.3%) and BioFreedom DCS (20.7%) for cardiac death (CD), myocardial infarction (MI), and stent thrombosis (ST)
- Significantly lower cd-TVR† with Resolute Onyx DES (4.8%) compared to BioFreedom DCS (7.1%) ($p = 0.035$)

*Third-party brands are trademarks of their respective owners.

†Endpoint were not powered.


Source: Windecker S, et al. Final Two-Year Results from the Randomized Onyx ONE Trial in High Bleeding Risk Patients Treated with 1-month DAPT. Presented at ACC 2021.



ONYX ONE CLEAR ANALYSIS

Reinforced Resolute Onyx™ DES was safe and effective in HBR patients on 1-month DAPT

BEAT PERFORMANCE GOAL DERIVED FROM CONTEMPORARY 1-MONTH DAPT TRIALS*



- Prospective, multicentre, single-arm analysis
- ~1500 patients included in primary endpoint analysis
- 1.6 high bleeding risk criteria per patient
- 4% BARC 3-5 bleeding rate at 1 year
- Primary endpoint results showed 7.0% cardiac death or myocardial infarction at one year, beating the performance goal of 9.7%
- Resolute Onyx was the first DES FDA indicated for HBR patients with 1-month DAPT labeling based on the results from this analysis

*ZEUS, LEADERS FREE, and SENIOR trials.
Source: Kandzari DE, et al. *Circ Cardiovasc Interv.* 2020;13:e009565.



Read the Onyx ONE Clear analysis article printed in *Circulation: Cardiovascular Interventions*




Medtronic

ACUTE CORONARY SYNDROME (ACS) SUBANALYSIS RESULTS

Reinforced Resolute Onyx™ DES was safe and effective in HBR, ACS patients on 1-month DAPT†

SIGNIFICANTLY LOWER MI IN HIGH ISCHEMIC RISK PATIENTS**



- Prespecified subanalysis from Onyx ONE global trial
- Compared Resolute Onyx DES to BioFreedom DCS
- 53% HBR ACS patients (n = 944)
- 1.7 average high bleeding risk criteria per patient
- Significantly higher device success with Resolute Onyx DES (92%) vs. BioFreedom DCS (87%) $p = 0.001$

*Third-party brands are trademarks of their respective owners.

†Results not adjusted for multiple comparisons.

**Endpoints were not powered.


Source: Kedhi A, et al. Outcomes in High Bleeding Risk Patients with Acute Coronary Syndrome with 1-Month DAPT: Insights from the Onyx ONE Trial. Presented at ESC Congress 2020.



ATRIAL FIBRILLATION (AF) SUBANALYSIS

Reinforced Resolute Onyx™ DES was safe and effective in HBR, AF patients on 1-month DAPT*

NO DIFFERENCE IN ISCHEMIC EVENTS AND BLEEDING (BARC 3–5)†



- Subanalysis from Onyx ONE Clear analysis
- AF patients are often on oral anticoagulants (OACs). Therefore, bleeding risk in patients on triple therapy (DAPT and OACs) is magnified.
- 36% HBR AF patients (n = 1491)
- Significantly higher average high bleeding risk criteria per AF patient (1.7) vs. no-AF patient (1.5, $p < 0.001$) with 87% of AF patients on oral anticoagulants

*Based on post-hoc analysis. Results not adjusted for multiple comparisons.

†Endpoint were not powered.

** P -values adjusted by propensity scores for baseline differences; all other p -values are not significant.

Source: Pasupati S, et al. Ischemic and Bleeding Outcomes in Patients With vs. Without Atrial Fibrillation: Analysis From the Onyx ONE Program. Presented at TCT Congress 2020.




Medtronic

COMPLEX PCI SUBANALYSIS

Reinforced Resolute Onyx™ DES was safe and effective in HBR, complex PCI patients on 1-month DAPT*

NO DIFFERENCE IN SAFETY AND EFFICACY DESPITE INCREASED LESION COMPLEXITY†



- Subanalysis from Onyx ONE Clear analysis
- Complex PCI patients met at least one characteristic¹:
 - 3 vessels treated
 - ≥ 3 lesions treated
 - Total stent length > 60 mm
 - Bifurcation with ≥ 2 stents implanted
 - Use of any atherectomy device
 - Left main surgical bypass graft[†]
 - Chronic total occlusion
- 1.6 high bleeding risk criteria per patient
- Complex PCI group included significantly higher multivessel disease (78.3% vs. 39.5% $p < 0.001$) and B2/C lesions (84.2% vs. 75.5% $p < 0.001$)

*Based on post-hoc analysis. Results not adjusted for multiple comparisons.

†Endpoint. Not powered.

**Propensity score adjusted p-values for differences in baseline characteristics.

Source: Kandzari D, et al. Complex PCI with 1-month DAPT in High Bleeding Risk Patients: Analysis from the Onyx ONE Clear Study. Presented at TCT 2020.


¹Dangas G, et al. *J Am Coll Cardiol*. 2020;75:2414-2424.



SEX SUBANALYSIS

Reinforced Resolute Onyx™ DES was safe and effective in HBR, female or male patients on 1-month DAPT*

LOW EVENT RATES INCLUDING CD/MI AND ST WITH NO SEX-BASED DIFFERENCES OTHER THAN TVR†



- Prespecified subanalysis from Onyx ONE Clear analysis
- 32% female population
- Significantly higher bleeding risk for females (1.6 average) vs. males (1.5 average) $p = 0.02$
- Differences in patient history and CAD symptoms highlight the need for sex-based analyses: ~28% non-STEMI in female patients but significantly higher prior interventions in men (previous PCI, CABG)

*Results not adjusted for multiple comparisons.

†Endpoint were not powered.

**Propensity score adjusted p-values for differences in baseline characteristics.


Source: Mehran, R et al. Sex-based Outcomes after PCI in Complex High Bleeding Risk Patients: Results from the Onyx ONE Clear Trial. Presented at SCAI 2021.



DIABETIC SUBANALYSIS

Reinforced Resolute Onyx™ DES was safe and effective in HBR patients with diabetes on 1-month DAPT*

NO DIFFERENCE IN SAFETY AND EFFICACY OBSERVED IN DM VS. NO-DM GROUPS DESPITE INCREASED LESION COMPLEXITY†



- Prespecified subanalysis from Onyx ONE Clear analysis
- Patients with diabetes mellitus (DM) historically show a high risk for ischemic events
- 40% DM patients
- No difference in average high bleeding risk criteria per patient between DM (1.6) and no-DM (1.5, $p = 0.06$)
- DM patients experience higher comorbidities prior to interventions, and longer stent lengths

*Results not adjusted for multiple comparisons.

†Endpoint were not powered.

**Propensity score adjusted p -values for differences in baseline characteristics.

Source: Kedhi E, et al. Diabetic High Bleeding Risk Patients with One-Month DAPT: Onyx ONE Clear Results. Presented at Euro PCR 2021.



Medtronic

SHORT-DAPT TRIALS COMPARISON

TRIAL DESIGN					LESION CHARACTERISTICS			PATIENT COMPLEXITY	
Name	Size	Comparator	Randomised Control Trial	Angiographic Exclusions	B2/C	Average Stented Length	Moderate/ Severe Calcified Lesions	ACS	Prior MI
ONYX ONE GLOBAL TRIAL¹	1996	BioFreedom™ DCS	✓	None	80%	38 mm	46%	53%	26%
ONYX ONE CLEAR ANALYSIS²	1506	Objective Performance Criteria (OPC)	✗	None	79%	37 mm	50%	49%	26.3%
XIENCE 28 STUDY³	1392	Single-arm Historical Control (XIENCE V USA Study 2008–2011)	✗	Excluded: Left main CTO ISR Overlapping stents SVG	36%	27 mm	Not Reported	34% No STEMI	15.8%
POEM (Synergy™ Stent)⁴	443	Objective Performance Criteria	✗	None	49%	24 mm	Not Reported	41%	Not Reported

*Third-party brands are trademarks of their respective owners.

¹Results include the Resolute Onyx DES treated arm.

²Windecker S, et al. *N Engl J Med*. 2020;382:1208-1218.

³Kandzari D, et al. *Circ Cardiovasc Interv*. 2020;13: e009565.



⁴Mehran R, et al. The Xience Short DAPT Program: Xience 90/28. Evaluating the Safety of 3-month and 1-month DAPT in HBR Patients. Presented at TCT Congress 2020.

⁵Stefanini G, et al. The POEM Study: One-Month DAPT in HBR Patients. Presented at PCR 2021.



WHY RESOLUTE ONYX™ DES FOR HBR PATIENTS ON 1-MONTH DAPT?

RESOLUTE ONYX DES IS DIFFERENT BY DESIGN TO PROMOTE FAST HEALING



Single-wire design

provides a fluid range of motion and the conformability needed for superior strut apposition¹

BioLinx™ biocompatible polymer

provides superior thromboresistance²

Fast healing occurs as evidenced by nearly 90% strut coverage at 30 days,³ **providing the option to shorten DAPT**

*Third-party brands are trademarks of their respective owners.
¹ Data on file at Medtronic.
² Jinnouchi H, et al. *Int J Cardiol.* 2021;327:52-57.
³ Roleder T, et al. *Postepy Kardiol Interwencyjnej.* 2019;15:143-150.



CE
2797

Europe
Medtronic Intl. Trading SARL
Tel: 41.21.802.7000

Asia Pacific
Medtronic Intl. Ltd.
Tel: 65.6436.5000

Latin America
Medtronic USA, Inc.
Tel: 786.709.4200

[medtronic.com/OnyxOne](https://www.medtronic.com/OnyxOne)

UC202200237 ML ©2021 Medtronic. All rights reserved. Medtronic and the Medtronic logo are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. For distribution only in markets where the Resolute Onyx coronary stent has been approved. Not for distribution in Japan, France, Canada, or Australia. 07/2021



Medtronic